



# PLC-Heating Chiller Manual

## I.Main screen



**Log in:** Click login in the upper right corner to pop up the user login interface. Only after the user logs in can you operate and modify the parameters. The system sets 8 users.

User name: 1, initial password: 111;

User name: 2, initial password: 222;

User name: 3, initial password: 333;

User name: 4, initial password: 444;

User name: 5, initial password: 555;

User name: 6, initial password: 666;

User name: 7, initial password: 777;

User name: 8, initial password: 12121;

The user name will be displayed if the user logs in successfully;

User 8 has the highest authority and can set all parameters;

Users 1-7 can set other parameters other than "Export High Temperature Protection",

"Export Low Temperature Protection", "Export High Temperature Protection",

"Heating P", "Heating I", and "Heating D";

You cannot enter the parameter setting interface without logging in.

**Quit:** Click "Quit" to exit interface

**Running time:** After the program mode is started, the user's running time will be displayed; each time it is started, the time will be restarted

**Current stage:** Display the current stage of the program mode

**Keep warm time:** Display the time that the current stage of the program mode has entered the heat preservation stage. Unit: minute

**Remaining holding time:** Display the remaining holding time of the current stage of the program mode. Unit: minute

**Heating output:** Display heating percentage, full power is 100%

**Outlet temperature:** The real-time temperature of the equipment outlet

**Material temperature:** The real-time temperature of the material in the reactor

**Target temperature:** The target temperature of the current control medium. If "Outlet temperature control" is selected, the set target temperature is the target temperature of the outlet;

If "material temperature control" is selected, the set target temperature is the target temperature of the material;

If it is in the program mode and the temperature rise time is set, the target temperature will automatically change once per minute

**Reverse phase:** white indicator light is normal, blue indicator light is phase loss or reverse phase failure

**Valve:** The white indicator light is normal, the blue indicator light is the valve has started

**Constant temperature:** The white indicator light is normal, the blue indicator light has entered the constant temperature, and the target temperature difference of  $\pm 0.2$  degrees will enter the constant temperature

**Circulation:** The white indicator light is normal, the blue indicator light is the circulating pump has started

**Heating:** White indicator light is normal, blue indicator light is heating has started

**Refrigeration:** The white indicator light is normal, the blue indicator light is the compressor has started


**Program mode:** This button is the automatic start button for stage heating. After starting, it will run according to the heating and holding time set by the user. When the program runs the current stage constant temperature time and detects that the next stage constant temperature time is 0, the program automatically stops

**Stand-alone mode:** This button does not execute the program to raise and lower the temperature, and it will keep the target temperature set by the user.

Outlet temperature control: select this button to set the outlet temperature as the target temperature

**Material temperature control:** Choose this button to set the material temperature as the target temperature

## II. Parameter setting



1. **Oil temperature correction:** when the displayed oil temperature does not match the actual temperature, it can be corrected to the same

2. **Material correction:** when the displayed material temperature does not match the actual temperature, it can be corrected to the same

3. **Start refrigeration temperature difference:** when the actual temperature is higher than the target temperature, the system will automatically start the compressor refrigeration

4. **Stop cooling temperature difference:** when the actual temperature is lower than the target temperature, the system will automatically stop the compressor cooling

5. **Oil temperature difference:** When the system takes the material temperature as the control target, in order to reach the target temperature faster, the oil temperature will be lower or higher than the temperature difference value of the material temperature, but

when the oil temperature is higher or lower than the target temperature in the kettle. When the temperature difference is set, it will stop heating or cooling until it returns to the temperature difference range to continue.

**6.Solenoid valve temperature:** When the oil temperature is higher than this value, the solenoid valve will start, the initial value is 40 degrees, do not change at will.

**7.Heating delay:** In order to reduce the system shock, the compressor will stop and the heating will not start within this time range.

**8.Temperature coefficient in the tank:** the time it takes for the materials in the kettle to rise or drop by one degree when the oil temperature reaches the highest value, unit s.

**9.Early action coefficient:** This coefficient divided by the temperature coefficient in the kettle is the temperature for heating or cooling in advance.

**10.Heating constant temperature difference:** When the material is selected for temperature control and heating, the oil temperature in the constant temperature stage is greater than the value of the target temperature of the material.

**11.Refrigeration constant temperature difference:** When the material is selected for temperature control and is cooling, the oil temperature in the constant temperature stage is lower than the value of the target temperature of the material.

**12.Delay protection:** protect the compressor, prevent the compressor from starting frequently, the factory setting is 180s.

**14.Export high temperature protection:** the highest value of oil temperature.

**15.Export low temperature protection:** the lowest value of oil temperature.

**16.Constant temperature alarm temperature difference:** when the temperature difference

in the constant temperature stage exceeds this value, it will alarm

**17.Heating P:** L factory value is 59.8, H factory value is 15.0, L is the low temperature section data, H is the high temperature section data

**18.Heating I:** L factory value is 3.67, H factory value is 3.67, L is the low temperature section data, H is the high temperature section data

**19.Heating D:** L factory value is 0, H factory value is 0, L is the low temperature section data, H is the high temperature section data


**20.Restore factory settings:** Click this button to restore all parameters to factory values

**21.System time calibration:** When the system time does not match the current time, you can click to calibrate, and the format must be the same, otherwise the calibration will fail

**22.Chinese:** switch the system display language to Chinese


**23.English:** The system display language is switched to English

### III. Trend curve



This interface displays the temperature curve and data of the last 8 hours, you can view the temperature value at a specific point in time, the up and down arrows can switch between the oil temperature and the material temperature

## IV. Manual screen



**Manual mode:** When starting the heating or compressor in the secondary mode, you need to start the circulating pump first

**Free mode:** When starting heating or compressor in secondary mode, there is no need to start the circulating pump

**Circulating pump:** Press this button to start the circulating pump, click again to stop the circulating pump

**Heating:** Press this button to start heating, click again to stop heating

**Compressor:** Press this button to start the compressor, click again to stop the compressor

**Manual trip:** Press this button, the trip will start and cut off the system power

# V.Stage setting

**SIEMENS**
SMART LINE

8/9/2021 5:01:36 PM
8

Running time: 0 D 0 H 0 M 0 S
Log in Quit

Main Screen

Setting

Trend Curve

Manual Screen

Phase Setting

Alarm

Phase	Heating/Cooling Time(min)	Insulation time(min)	Target Tem (°C)	
Phase1	<input type="text" value="0"/>	<input type="text" value="0"/>	<input type="text" value="0.0"/>	
Phase2	<input type="text" value="0"/>	<input type="text" value="0"/>	<input type="text" value="0.0"/>	Phase1-10
Phase3	<input type="text" value="0"/>	<input type="text" value="0"/>	<input type="text" value="0.0"/>	
Phase4	<input type="text" value="0"/>	<input type="text" value="0"/>	<input type="text" value="0.0"/>	Phase11-20
Phase5	<input type="text" value="0"/>	<input type="text" value="0"/>	<input type="text" value="0.0"/>	
Phase6	<input type="text" value="0"/>	<input type="text" value="0"/>	<input type="text" value="0.0"/>	Phase21-30
Phase7	<input type="text" value="0"/>	<input type="text" value="0"/>	<input type="text" value="0.0"/>	
Phase8	<input type="text" value="0"/>	<input type="text" value="0"/>	<input type="text" value="0.0"/>	
Phase9	<input type="text" value="0"/>	<input type="text" value="0"/>	<input type="text" value="0.0"/>	Phase31-40
Phase10	<input type="text" value="0"/>	<input type="text" value="0"/>	<input type="text" value="0.0"/>	

Save
Previous Group
0
Next Group
Apply

TOUCH

**SIEMENS**
SMART LINE

8/9/2021 5:01:50 PM
8

Running time: 0 D 0 H 0 M 0 S
Log in Quit

Main Screen

Setting

Trend Curve

Manual Screen



Phase Setting

Alarm

Phase	Heating/Cooling Time(min)	Insulation time(min)	Target Tem (°C)	
Phase11	<input type="text" value="0"/>	<input type="text" value="0"/>	<input type="text" value="0.0"/>	
Phase12	<input type="text" value="0"/>	<input type="text" value="0"/>	<input type="text" value="0.0"/>	Phase1-10
Phase13	<input type="text" value="0"/>	<input type="text" value="0"/>	<input type="text" value="0.0"/>	
Phase14	<input type="text" value="0"/>	<input type="text" value="0"/>	<input type="text" value="0.0"/>	Phase11-20
Phase15	<input type="text" value="0"/>	<input type="text" value="0"/>	<input type="text" value="0.0"/>	
Phase16	<input type="text" value="0"/>	<input type="text" value="0"/>	<input type="text" value="0.0"/>	Phase21-30
Phase17	<input type="text" value="0"/>	<input type="text" value="0"/>	<input type="text" value="0.0"/>	
Phase18	<input type="text" value="0"/>	<input type="text" value="0"/>	<input type="text" value="0.0"/>	
Phase19	<input type="text" value="0"/>	<input type="text" value="0"/>	<input type="text" value="0.0"/>	Phase31-40
Phase20	<input type="text" value="0"/>	<input type="text" value="0"/>	<input type="text" value="0.0"/>	

Save
Previous Group
0
Next Group
Apply

TOUCH



The system can set 40 sections of programs to raise and lower the temperature, and 6

groups of formulas can be set, which is convenient for users to call.

**Temperature rise and fall time:** This parameter sets the heating or cooling time of each stage. The system will change the target temperature every minute according to the current temperature difference. After the temperature rise and fall time, the system will move to the final target temperature of the stage. When the constant temperature time goes After the system has not reached the target temperature, the system will continue to the target temperature at this time. When it reaches the target temperature by 0.2 degrees, it will enter the constant temperature timer; when the temperature rise and fall time is not required, it can be set to 0, and the system directly rises and falls to the final target temperature.

**Constant temperature time:** This parameter sets the constant temperature time of each stage, and the constant temperature timer starts when the current temperature reaches the target temperature  $\pm 0.2$  degrees. After the current stage constant temperature time is over, continue to the next stage, when the system judges that the next stage constant temperature time is 0, it will automatically stop.

**Save:** After setting, click the save button to save the set parameters

**Previous Group:** Click this button to call up the last set of formula parameters




: You can also directly call up the corresponding formula group by directly inputting the value

**Next group:** Click this button to call up the next set of formula parameters

**0:** Click this button to clear the current formula parameters

## VI.Alarm information



This interface will display current and historical alarm information, click 0 to clear the historical alarm information